

Influenza A virus matrix protein antibody [B313M]

Cat. No. GTX41302

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, EIA, Lateral Flow
Reactivity	Influenza A virus

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
EIA	Assay dependent
Lateral Flow	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	4.49 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Influenza A Virus A/Puerto Rico/8/34 (H1N1) and A/Bangkok/1/79 (H3N2) virus.
Purification	Protein A purified From tissue culture supernatant
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).