

Infectious Bursal Disease virus (IBDV) VP2 antibody [B826M]

Cat. No. GTX41332

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA, Sandwich ELISA
Reactivity	Infectious Bursal Disease virus

Package

500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Not tested in other applications.

Product Note Detects VP2 protein of IBD virus. Does not crossreact with other avian viruses.

Properties

Form	Liquid
Buffer	PBS pH7.4
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	6.20 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified virulent IBD virus (field strain) derived from infected chicken bursas.
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2025 / 05 / 02 Page 1 of 1