

Dehydroepiandrosterone antibody [DHEA(-3):1-19.1 (BGN/117/1191)]

Cat. No. GTX41444

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA
Reactivity	Species independent

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	1/1000-1/10,000

Not tested in other applications.

Product Note This antibody detects dehydroepiandrosterone (DHEA) a steroid prohormone which is the precursor of all sex steroids.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	DHEA-3-HS-BSA conjugate
Purification	Protein A purified From tissue culture supernatant
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX41444 ELISA Image

ELISA analysis of Dehydroepiandrosterone using GTX41444 Dehydroepiandrosterone antibody [DHEA(-3):1-19.1 (BGN/117/1191)].



For full product information, images and publications, please visit our [website](https://www.genetex.com).