

Dectin-2 antibody [D2.11E4]

Cat. No. GTX41453

Host	Rat
Clonality	Monoclonal
Isotype	lgG2a
Applications	IHC-Fr, FCM, IP
Reactivity	Mouse

References (1)
Package
25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-Fr	Assay dependent
FCM	1/25-1/50
IP	Assay dependent

Note : Use 10 μ l of the suggested working dilution to label 10 6 cells in 100 μ l.

Not tested in other applications.

Properties	
Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Dectin-2 (extracellular domain) Fc Fusion protein.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 03 Page 1 of 2



DATA IMAGES



GTX41453 IHC-Fr Image

IHC analysis of mouse spleen tissue using GTX41453 Dectin-2 antibody [D2.11E4].



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 03 Page 2 of 2