

Cryptosporidium antibody [BEL 0126]

Cat. No. GTX41587

Host	Mouse
Clonality	Monoclonal
Isotype	IgG3
Applications	ICC/IF, FCM, ELISA
Reactivity	Cryptosporidium

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1/50-1/200
FCM	Assay dependent
ELISA	1/2,000-1/8,000

Not tested in other applications.

Product Note This antibody recognizes a membrane antigen expressed by the oocysts of Cryptosporidium sp.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified cryptosporidium oocysts from bovine faeces.
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 1