

Hepatitis C virus Core Antigen antibody [11-B3]

Cat. No. GTX41722

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-Fr, ELISA
Reactivity	Hepatitis C virus

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1,000
ICC/IF	Assay dependent
IHC-Fr	1:20
ELISA	1:1,000,000

Not tested in other applications.

Product Note HCV Core Antigen (Amino Acids 70-90).

Properties


Form	Liquid
Buffer	PBS
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant HCV core antigen (genotype 1b).
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note
 Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX41722 ICC/IF Image

The data was published in the journal PLoS One in 2010. [PMID: 20877724](https://doi.org/10.1371/journal.pone.0208772)



For full product information, images and publications, please visit our [website](https://www.genetex.com).