

## Hepatitis B virus Surface Antigen antibody [B521M]

Cat. No. GTX41736

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, Lateral Flow, Sandwich ELISA
Reactivity	Hepatitis B virus

References ( 1 )

Package

1 mg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Lateral Flow	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture : GTX41736, detection : GTX41745**

Not tested in other applications.

**Product Note** Hepatitis B Surface Antigen, ad and ay subtypes.

## Properties

Form	Liquid
Buffer	Borate-buffered saline
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	8.67 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified Hepatitis B Surface Antigen.
Purification	Purified by Ion Exchange Chromatography
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).



For full product information, images and publications, please visit our [website](#).