

## CD70 antibody [Bu69]

Cat. No. GTX41860

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-Fr, FCM
Reactivity	Human

## Package

200 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-Fr	Assay dependent
FCM	1/50-1/200

**Note : Use 10µl of the suggested working dilution to label 10<sup>6</sup> cells in 100µl.**

Not tested in other applications.

## Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX41860 IHC-Fr Image**


IHC-Fr analysis of a human tonsil with Mouse anti Human CD70 antibody, clone Bu69 (GTX41860), red in A and Mouse anti Human CD21 antibody, clone LB21 (GTX74728), green in B. C is the merged image with nuclei counterstained blue using DAPI. High power


**GTX41860 IHC-Fr Image**


IHC-Fr analysis of human tonsil tissue using GTX41860 CD70 antibody [Bu69].


**GTX41860 IHC-Fr Image**

IHC-Fr analysis of human tonsil tissue using GTX41860 CD70 antibody [Bu69].


**GTX41860 IHC-Fr Image**

IHC-Fr analysis of a human tonsil with Mouse anti Human CD70 antibody, clone Bu69 (GTX41860), red in A and Mouse anti Human CD21 antibody, clone LB21 (GTX74728), green in B. C is the merged image with nuclei counterstained blue using DAPI. Low power


**GTX41860 IHC-Fr Image**

IHC-Fr analysis of human tonsil tissue using GTX41860 CD70 antibody [Bu69].



For full product information, images and publications, please visit our [website](https://www.genetex.com).