

CD63 antibody [AD1]

Cat. No. GTX41877

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, FCM, IP, ELISA
Reactivity	Human, Mouse, Rat

References (7)

Package

25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
FCM	1/10-1/200
IP	Assay dependent
ELISA	Assay dependent

Note : Use 10µl of the suggested working dilution to label 10⁶ cells in 100µl.

Not tested in other applications.

Calculated MW 26 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	RBL-2H3 rat basophilic leukaemia cell line.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX41877 WB Image

The data was published in the journal Stem Cell Res Ther in 2019. [PMID: 30876481](#)



GTX41877 WB Image

The data was published in the journal Nat Commun in 2017. [PMID: 29038531](#)



GTX41877 FCM Image

The data was published in the journal Stem Cell Res Ther in 2019. [PMID: 30876481](#)



For full product information, images and publications, please visit our [website](#).