

ENaC Beta antibody

Cat. No. GTX41971

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IHC
Reactivity	Human, Mouse, Rat

References (1)
Package
100 µl

Applications

Application Note

We recommend a working dilution of 1:1000 (WB).

Calculated MW 72 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Rabbit immunoglobulin, 50% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	AA617-638 of the rat sequence
Purification	Affinity purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx41971 WB Image

Western blot analysis of ENaC beta in rat kidney tissues using a 1:1000 dilution of GTx41971.



GTx41971 ICC/IF Image

IF analysis of ENaC beta in rat kidney tissues using a dilution of GTx41971.



GTx41971 IHC-P Image

The data was published in the journal Front Physiol in 2018. [PMID: 30050465](https://pubmed.ncbi.nlm.nih.gov/30050465/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).