

Nitrotyrosine antibody [39B6]

Cat. No. GTX41979

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF, IHC-P, IHC-Fr, IP, ELISA
Reactivity	Species independent

References (8)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1400
ICC/IF	Assay dependent
IHC-P	1:100
IHC-Fr	Assay dependent
IP	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Product Note Recognizes 3-nitrotyrosine moieties. No detectable cross-reactivity with non-nitrated tyrosine. Not species specific.

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	3-(4-hydroxy-3-nitrophenylacetamido) propionic acid-bovine serum albumin
Purification	Protein G purified
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx41979 WB Image

Nitrotyrosine visualized on a retinal injury model using GTx41979.


GTx41979 ICC/IF Image

IHC-P analysis of rat liver tissue using GTx41979 Nitrotyrosine antibody [39B6].


GTx41979 IHC-P Image

Nitrotyrosine tested on Bouin's fixed paraffin-embedded backskin sections of transgenic mice using a dilution of GTx41979



For full product information, images and publications, please visit our [website](#).