

TrpC5 antibody [S67-15]

Cat. No. GTX42023

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, IP, IHC
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

We recommend a working dilution of iug/mL (WB), 0.1ug/mL (IHC/ICC), 1-10ug/mL (IF).

Calculated MW 111 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide amino acids 827-845 of human TrpC5 (also known as short transient receptor potential channel 5, and Htrp5
Purification	Protein G Purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX42023 IHC Image

IHC analysis of TrpC5 tested on mouse brain using a dilution of GTX42023. Nuclear membrane staining of every cells.



For full product information, images and publications, please visit our [website](#).