

## CD252 antibody [ATM-2] (Low endotoxin, azide free)

Cat. No. GTX42204

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	FACS, Functional Assay
Reactivity	Rat

**Package**  
500 µg

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FACS	Assay dependent
Functional Assay	Assay dependent

Not tested in other applications.

**Calculated MW** 22 kDa. ( [Note](#) )

## PROPERTIES

Form	Liquid
Buffer	PBS
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	ACI-S1 strain rat T-cell lymphoma cell line.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).