

# Hepatitis C virus Core Antigen antibody [B056M]

**Cat. No. GTX42537**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Application</b>	ICC/IF, ELISA, Sandwich ELISA
<b>Reactivity</b>	Hepatitis C virus

**Package**  
500 µg

## APPLICATION

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Not tested in other applications.

**Product Note** Hepatitis C Virus (HCV) Core Antigen. Reactive with genotypes 1a and 2.

## PROPERTIES

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	0.10 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant HCV core antigen, genotype 1a.
<b>Purification</b>	Protein A purified
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).