

hCG beta antibody [BDI585]

Cat. No. GTX42540

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	Lateral Flow
Reactivity	Human

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

Lateral Flow	Assay dependent
--------------	-----------------

Not tested in other applications.

Product Note

Cross-Reactivity: (Tested as an antibody pair with a MAb to hCG alpha in a sandwich immunoassay to detect intact hCG.) beta hCG 100% Intact hCG 115% LH <1% FSH <1% TSH <1%

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	7.80 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified hCG beta
Purification	Purified by DEAE Chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).