

hCG alpha antibody [B181M]

Cat. No. GTX42569

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	ELISA
Reactivity	Human

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

Product Note Reacts with the alpha chain of FSH, LH, TSH, and hCG.

Properties

Form	Liquid
Buffer	10mM Tris, 150mM NaCl
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	3 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified hCG alpha chain.
Purification	Purified by Ion Exchange Chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).