

Foot and Mouth Disease virus, Serotype 01 antibody [B473M]

Cat. No. GTX42730

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ELISA, Immunodiffusion
Reactivity	Foot and Mouth Disease virus

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Immunodiffusion	Assay dependent

Not tested in other applications.

Product Note Specific for FMDV virulent serotype O1.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	2.10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Inactivated FMD Virus
Purification	Protein G purified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).