

Adiponectin antibody [Adn27]

Cat. No. GTX42852

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ELISA
Reactivity	Human

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human adiponectin purified from E.coli.
Purification	Protein A purified From tissue culture supernatant
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX42852 ELISA Image

Sandwich ELISA analysis of human adiponectin binding using Mouse anti Human adiponectin (GTX42852) as a capture reagent and biotinylated Mouse anti Human adiponectin (GTX49335) as a detection reagent. A plasma sample diluted 1:4 (green) is shown



For full product information, images and publications, please visit our [website](https://www.genetex.com).