

Goat Anti-Human lambda light chain antibody, F(ab')₂ fragment, pre-adsorbed(PE)

Cat. No. GTX42962

Host	Goat
Clonality	Polyclonal
Isotype	IgG F(ab') ₂
Applications	FCM
Reactivity	Human

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

FCM	1/20-1/30
-----	-----------

Note : Use 10µl of the suggested working dilution to label 10⁶ cells in 100µl.

Not tested in other applications.

Product Note

This antibody recognises human lambda light chains and has been cross absorbed against pooled human myeloma proteins with kappa light chains.

Properties

Form	Liquid
Buffer	PBS, Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.25 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified by affinity chromatography
Conjugation	Phycoerythrin (PE) Wavelength

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).