

## EBV Early Antigen-R Homologue Bcl-2 antibody [BDI502]

## Cat. No. GTX43014

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	WB, IHC-P, IP, ELISA
Reactivity	Epstein-Barr virus

Package 500 μg

## Applications

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IP	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

**Product Note** Recognizes Epstein-Barr virus. Specific for the early antigen (restricted) homologue Bcl-2.

Properties	
Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Infected Cell Lysate
Purification	Protein A purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 07 Page 1 of 1

€ 886-3-6208988 📻 886-3-6208989 🐷 infoasia@genetex.com