

# CD45R antibody [20.96] (PE)

# Cat. No. GTX43071

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	FCM
Reactivity	Sheep

Package 100 test

# Applications

#### **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Neat

Note : Use  $10\mu l$  of the suggested working dilution to label  $10^6$  cells in  $100\mu l$  .

Not tested in other applications.


Properties	
Form	Liquid
Buffer	PBS, 1% BSA, 5% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Ovine thymocytes.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Phycoerythrin (PE) Wavelength
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 10 / 26 Page 1 of 2

## DATA IMAGES



## GTX43071 FCM Image

FACS analysis of sheep peripheral blood lymphocytes using GTX43071 CD45R antibody [20.96] (PE).



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 10 / 26 Page 2 of 2