

# IL2 Receptor alpha antibody [OX-39] (PE)

**Cat. No. GTX43135**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Mouse      |
| <b>Clonality</b>    | Monoclonal |
| <b>Isotype</b>      | IgG1       |
| <b>Applications</b> | FCM        |
| <b>Reactivity</b>   | Rat        |

**Package**  
100 test

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| FCM                | Neat                 |

**Note : Use 10µl of the suggested working dilution to label 10<sup>6</sup> cells in 100µl.**

Not tested in other applications.

## Properties

|                      |                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                     |
| <b>Buffer</b>        | PBS, 1% BSA, 5% Sucrose                                                                                                    |
| <b>Preservative</b>  | 0.09% Sodium azide                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light. |
| <b>Concentration</b> | Batch dependent (Please refer to the vial label for the specific concentration.)                                           |
| <b>Immunogen</b>     | Stimulated Rat T cells                                                                                                     |
| <b>Purification</b>  | Protein G purified<br>From tissue culture supernatant                                                                      |
| <b>Conjugation</b>   | Phycoerythrin (PE) <a href="#">Wavelength</a>                                                                              |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX43135 FCM Image

FACS analysis of Con A activated rat splenic lymphocytes using GTX43135 IL2 Receptor alpha antibody [OX-39] (PE).



For full product information, images and publications, please visit our [website](#).