

Goat Anti-Human lambda light chain antibody, F(ab')₂ fragment, pre-adsorbed(HRP)

Cat. No. GTX43288

Host	Goat
Clonality	Polyclonal
Isotype	IgG F(ab') ₂
Application	WB, IHC-P, IHC-Fr, ELISA
Reactivity	Human

Package

1 ml

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
ELISA	1/4,000-1/8,000

Not tested in other applications.

Product Note

GTX43288 recognises human lambda light chains. This antibody has been cross absorbed against pooled human myeloma proteins with kappa light chains.

PROPERTIES

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Purification	Purified by affinity chromatography
Conjugation	Horseradish peroxidase(HRP)



For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).