

D-Dimer antibody [DD41] (Biotin)

Cat. No. GTX43437

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ELISA
Reactivity	Human

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

Product Note

Recognizes high molecular weight fibrin degradation products. Recognizes D-dimer in Western blot under non-reducing conditions. Interacts with beta-chain of D-dimer in Western blot under reducing conditions. Does not cross-react with fibrinogen and D-monomer.

Properties

Form	Liquid
Buffer	This detergent is supplied as 25% in water
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mixture of D-dimer and high molecular weight fibrin degradation products
Purification	Protein A purified
Conjugation	Biotin

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).