

CD8 antibody [38.65] (FITC)

Cat. No. GTX43502

| | |
|---------------------|---------------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG2a |
| Applications | FCM |
| Reactivity | Goat, Sheep, Bovine |

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| FCM | Neat-1/10 |

Note : Use 10µl of the suggested working dilution to label 10⁶ cells in 100µl.

Not tested in other applications.

Properties

| | |
|----------------------|--|
| Form | Liquid |
| Buffer | PBS, 1% BSA |
| Preservative | 0.09% Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light. |
| Concentration | 0.1 mg/ml (Please refer to the vial label for the specific concentration.) |
| Immunogen | Ovine efferent lymphocytes. |
| Purification | Protein G purified From tissue culture supernatant |
| Conjugation | Fluorescein isothiocyanate (FITC) Wavelength |


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note
Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

**GTX43502 FCM Image**

FACS analysis of sheep peripheral blood lymphocytes using GTX43502 CD8 antibody [38.65] (FITC).



For full product information, images and publications, please visit our [website](#).