

Trident ELISA Coating Buffer (5X)

Cat. No. GTX436344
Applications ELISA

Package
50 ml

PRODUCT

Summary Trident ELISA Coating Buffer is used for immobilization of an antigen or antibody onto the microplate surface.

Applications

Application Note

The GeneTex Trident ELISA Coating buffer (5X) must be diluted 5-fold using double-distilled water (ddH₂O) before use.

Properties

Form Liquid

Preservative 0.05% ProClin 300

Storage Store at 2-8°C.

Note For *In vitro* laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES


GTX436344 Image


For full product information, images and publications, please visit our [website](https://www.genetex.com).