

# Trident ELISA Coating Buffer (5X)

**Cat. No. GTX436344**

**Applications** ELISA

**Package**  
50 ml

## PRODUCT

**Summary** Trident ELISA Coating Buffer is used for immobilization of an antigen or antibody onto the microplate surface.

## Applications

## Application Note

The GeneTex Trident ELISA Coating buffer (5X) must be diluted 5-fold using double-distilled water (ddH<sub>2</sub>O) before use.

## Properties

**Form** Liquid

**Preservative** 0.05% ProClin 300

**Storage** Store at 2-8°C.

**Note** For *In vitro* laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

## DATA IMAGES



**GTX436344 Image**



For full product information, images and publications, please visit our [website](https://www.genetex.com).