

## Trident ELISA Coating Buffer (5X)

Cat. No. GTX436344

Application	ELISA
-------------	-------

Package	50 ml
---------	-------

## PRODUCT

Summary	Trident ELISA Coating Buffer is used for immobilization of an antigen or antibody onto the microplate surface.
---------	--

## APPLICATION

## Application Note

The GeneTex Trident ELISA Coating buffer (5X) must be diluted 5-fold using double-distilled water (ddH<sub>2</sub>O) before use.

## PROPERTIES

Form	Liquid
------	--------

Storage	Store at 2-8°C.
---------	-----------------

Note	For <i>In vitro</i> laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	--

## DATA IMAGES



GTX436344 Image



For full product information, images and publications, please visit our [website](#).

Date 2024 / 05 / 09 Page 1 of 1