

Trident ELISA Blocking Buffer (5X)

Cat. No. GTX436345

Applications ELISA**Package**
50 ml

PRODUCT

Summary

Trident ELISA Blocking Buffer is used to prevent sample component binding directly to the microplate surface as well as to stabilize bound biomolecules and decrease nonspecific interactions.

Applications

Application Note

The GeneTex Trident ELISA Blocking buffer (5X) must be diluted 5-fold using double-distilled water (ddH₂O) before use.

Properties

Form Liquid**Preservative** 0.05% ProClin 300**Storage** Store at 2-8°C.**Note** For *In vitro* laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES



GTX436345 Image



For full product information, images and publications, please visit our [website](#).

Date 2026 / 02 / 01 Page 1 of 1