

Trident ELISA Diluent Buffer (5X)

Cat. No. GTX436346**Applications** ELISA**Package**
50 ml**PRODUCT****Summary**

Trident ELISA Diluent Buffer is suitable for diluting samples (e.g., purified antibodies, proteins, human serum/plasma, cell culture supernatant) and detection reagents that utilize peroxidase activity.

Applications**Application Note**

The GeneTex Trident ELISA Diluent buffer (5X) must be diluted 5-fold using double-distilled water (ddH₂O) before use.

Properties**Form** Liquid**Preservative** 0.05% ProClin 300**Storage** Store at 2-8°C.**Note** For *In vitro* laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.**DATA IMAGES****GTX436346 Image**

For full product information, images and publications, please visit our [website](#).