

Trident ELISA Diluent Buffer (5X)

Cat. No. GTX436346
Applications ELISA

Package
50 ml

PRODUCT

Summary

Trident ELISA Diluent Buffer is suitable for diluting samples (e.g., purified antibodies, proteins, human serum/plasma, cell culture supernatant) and detection reagents that utilize peroxidase activity.

Applications

Application Note

The GeneTex Trident ELISA Diluent buffer (5X) must be diluted 5-fold using double-distilled water (ddH₂O) before use.

Properties

Form Liquid

Preservative 0.05% ProClin 300

Storage Store at 2-8°C.

Note For *In vitro* laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES


GTX436346 Image


For full product information, images and publications, please visit our [website](https://www.genetex.com).