

Cystatin C antibody [Cyst-23]

Cat. No. GTX43872

| | |
|---------------------|-----------------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Applications | ELISA, Sandwich ELISA |
| Reactivity | Human |

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| ELISA | Assay dependent |
| Sandwich ELISA | Assay dependent |

Note : Capture : GTX43872, Detection : GTX43943 / GTX43944

Not tested in other applications.

Properties

| | |
|----------------------|---|
| Form | Liquid |
| Buffer | PBS |
| Preservative | 0.09% Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| Concentration | 6.50 mg/ml (Please refer to the vial label for the specific concentration.) |
| Immunogen | Cystatin C purified from human urine. |
| Purification | Protein A purified |
| Conjugation | Unconjugated |

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).