

C4 antibody

Cat. No. GTX44051

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	Immunoelectrophoresis, Turbidimetry
Reactivity	Human

Package
10 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Immunoelectrophoresis	Assay dependent
Turbidimetry	Assay dependent

Not tested in other applications.

Product Note Monospecific by Immunoelectrophoresis (IEP) against normal human plasma and 2X concentrated normal human serum.

Properties

Form	Liquid
Buffer	50mM Tris, 500mM NaCl
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified Human C4
Purification	Purified by Delipidation followed by solid phase absorption.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).