

Citrate Antigen Retrieval Solution (10X, pH 6.0) HistoMAX™

Cat. No. GTX440621

Applications IHC-P

Package
100 ml

Applications

Application Note

This buffer is at 10X concentration. It must be diluted 1:10 to obtain the 1X (10 mM Citrate) working buffer.

Properties

Form Liquid

Preservative No preservatives

Storage Store at 2-8°C, DO NOT FREEZE.


Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX440621 IHC-P Image



For full product information, images and publications, please visit our [website](https://www.genetex.com).

Date 2026 / 10 / 07 Page 1 of 1