

## Tris-EDTA Antigen Retrieval Solution (10X, pH 9.0) HistoMAX™

Cat. No. GTX440623

**Applications** IHC-P**Package**  
100 ml

## Applications

## Application Note

This buffer is at 10X concentration. It must be diluted 1:10 to obtain the 1X (10 mM Citrate) working buffer.

## Properties

**Form** Liquid**Preservative** No preservatives**Storage** Store at 2-8°C, DO NOT FREEZE.**Conjugation** Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



GTX440623 IHC-P Image



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 1