

Collagen I antibody

Cat. No. GTX44085

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, IHC-Fr, ELISA, RIA
Reactivity	Bovine

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:500
IHC-Fr	1:40
ELISA	1:200
RIA	Assay dependent

Not tested in other applications.

Product Note

Cross-reactivity (RIA):Bovine Collagen type I 100% Bovine Collagen types II, IX, XI <0.1% Bovine Fibronectin <0.1% Human Collagen Type I <2.0% Chicken Collagen Type I <0.1%

Properties

Form	Liquid
Buffer	saline buffer
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Collagen type I extracted and purified from bovine skin
Purification	Purified by Ion Exchange Chromatography (DEAE-Trisacryl)
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).