

Chlamydia trachomatis MOMP antibody [135]

Cat. No. GTX44138

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	Immunoassay
Reactivity	Chlamydia

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Immunoassay	Assay dependent
Not tested in other applications.	

Product Note Reactive with the following serovars: A, B, Ba, C, D, E, F, G, H, I, J, K, L1, L2, L3.

Properties

Form	Liquid
Buffer	Ascites
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	C. trachomatis elementary bodies.
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).