

## Respiratory Syncytial virus antibody [BDI065]

## Cat. No. GTX44267

Host	Mouse	
Clonality	Monoclonal	
Isotype	lgG2a	
Applications	ICC/IF, ELISA, Lateral Flow, Neutralizing/Inhibition, Sandwich ELISA	
Reactivity	Respiratory syncytial virus	

References ( 1 ) Package 500 µg

## Applications

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution	
ICC/IF	Assay dependent	
ELISA	Assay dependent	
Lateral Flow	Assay dependent	
Neutralizing/Inhibition	Assay dependent	
Sandwich ELISA	Assay dependent	
Note: Capture: GTX39261/GTX39258, Detection: GTX44267		

Not tested in other applications.

**Product Note** Specific for the fusion protein of RSV, types A & B. Neutralizes.

Properties	
Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	RSV viral lysate (Long strain).
Purification	Protein A purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 12 Page 1 of 2

€ 886-3-6208988 📻 886-3-6208989 🐷 infoasia@genetex.com



For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 12 Page 2 of 2