

Respiratory Syncytial virus antibody [BDI065]

Cat. No. GTX44267

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ICC/IF, ELISA, Lateral Flow, Neutralizing /Inhibition, Sandwich ELISA
Reactivity	Respiratory syncytial virus

References (1)

Package

500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
ELISA	Assay dependent
Lateral Flow	Assay dependent
Neutralizing /Inhibition	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture : GTX39261/GTX39258, Detection : GTX44267

Not tested in other applications.

Product Note Specific for the fusion protein of RSV, types A & B. Neutralizes.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	RSV viral lysate (Long strain).
Purification	Protein A purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).