

## Campylobacter jejuni antibody [BDI727]

## Cat. No. GTX44282

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	WB, ELISA, Lateral Flow
Reactivity	Campylobacter jejuni

Package 500 μg

## Applications

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Lateral Flow	Assay dependent
Not tosted in other applications	

Not tested in other applications.

Product Note

Specific to C. jejuni in whole cell ELISA. Reactive with a ~100 kDa antigen in Western Blot. Does notreact with H. pylori, E. coli (several), Salmonella (several), P. aeruginosa, Yersinia and Citrobacter.

Properties	
Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1.60 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	C. jejuni whole organism.
Purification	Protein A purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 08 Page 1 of 1