

CA125 antibody [X306]

Cat. No. GTX44294

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA, Sandwich ELISA
Reactivity	Human

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture : GTX44294, Detection : GTX10029/GTX10033/GTX21107

Not tested in other applications.

Calculated MW 1519 kDa. ([Note](#))

Product Note Epitope specificity group A (ISOBM classification) similar to OC125.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09 % Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	CA125 antigen purified from human ovarian carcinoma (MW > 1 MDa).
Purification	Purified by Protein A
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).