

BrdU antibody [131-14871(AH4H7-1)]

Cat. No. GTX44318

Host	Mouse	
Clonality	Monoclonal	
Isotype	lgG1	
Applications	FCM, IHC	
Reactivity	Species independent	

Package $100 \, \mu g$

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Bromodeoxyuridine (BrdU), a thymidine analogue, is incorporated into DNA during replication. Anti-BrdUreacts with cells **Product Note** that incorporated BrdU during proliferation activity in a variety of specimens such astumors, blood, cell lines and tissues.

Properties	
Form	Liquid
Buffer	10mM Phosphate, 150mM NaCl
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Protein A purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our website.

Date 2025 / 12 / 13 Page 1 of 1