

BNP antibody

Cat. No. GTX44337

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, ELISA, RIA
Reactivity	Human, Rat

Package
25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:1000
ELISA	Assay dependent
RIA	Assay dependent

Not tested in other applications.

Product Note Specific for human and rat Brain Natriuretic Peptide (BNP).

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic human Brain Natriuretic Polypeptide covalently attached to a carrier protein (bovine
Purification	Purified by Precipitation, chromatography and dialysis.
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).