

Biotin antibody [B169M]

Cat. No. GTX44344

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, EIA, Lateral Flow
Reactivity	Species independent

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
EIA	Assay dependent
Lateral Flow	Assay dependent

Not tested in other applications.

Product Note Recognizes free biotin, biotin conjugated to KLH or HRP, and biotinylated DNA probes.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	8.85 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Biotin-KLH
Purification	Purified by Ion Exchange Chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).