

Benzoylecgonine (Cocaine) antibody [B852M]

Cat. No. GTX44362

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1

Package
1 mg

PROPERTIES

Form	Liquid
Buffer	PBS pH7.0
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	4.30 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	BZE- BSA Conjugate
Purification	Purified IgG
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).