

Apolipoprotein A1 antibody

Cat. No. GTX44402

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	RIA
Reactivity	Mouse

Package
1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
RIA	Assay dependent

Not tested in other applications.

Product Note Recognizes mouse and rat Apo AI. Does not react well with human Apo AI.

Properties

Form	Liquid
Buffer	Serum, 100 U/ml Penicillin, 10µg/ml Streptomycin, 25 ng/ml Amphotericin B
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Mouse Apolipoprotein AI from pooled mouse plasma high density lipoprotein.
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).