

Adiponectin antibody [B871M]

Cat. No. GTX44469

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IP, ELISA, Sandwich ELISA
Reactivity	Human

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IP	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX44468, Detection: GTX44469 or Capture: GTX44469, Detection: GTX44468

Please notice that detection antibody needs to be conjugated to HRP to function when paired with capture antibody. Please contact us for custom HRP-conjugated antibody.

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	3.30 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Native Human Adiponectin
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).