

Adenovirus hexon antibody [B934M]

Cat. No. GTX44483

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, ELISA, Lateral Flow, Sandwich ELISA
Reactivity	Adenovirus

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
ELISA	Assay dependent
Lateral Flow	Assay dependent
Sandwich ELISA	Assay dependent

Not tested in other applications.

Product Note

Recognizes the hexon group antigen of Adeno types 1, 2, 3, 4, 5, 6, 7a, 8, 31, 40 and 41. Not reactive with Influenza A, Influenza B, RSV, Para 1, 2 and 3, Mycoplasma pneumoniae, H. pylori and mammalian cells.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	2.40 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Partially purified Adenovirus.
Purification	Protein A purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).