

Adenovirus hexon antibody [B315M]

Cat. No. GTX44484

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB, EIA, Lateral Flow
Reactivity	Adenovirus

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
EIA	Assay dependent
Lateral Flow	Assay dependent

Not tested in other applications.

Product Note Adenovirus hexon protein (genus specific epitope). Recognizes all human serotypes.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	6.68 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Adenovirus type 3 (ATCC strain VR847).
Purification	Protein A purified From tissue culture supernatant
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).