

ATP6V0A1 antibody, C-term

Cat. No. GTX44653

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Mouse, Pig

References (2)

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 96 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a C-terminal region of Mouse ATP6V0A1.
Purification	Affinity Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX44653 WB Image

WB analysis of mouse uterus tissue using GTX44653 Atp6v0a1 antibody at 0.2-1µg/ml.



GTX44653 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX44653 Atp6v0a1 antibody at 1:250.

A. untransfected HeLa cells

B. mATP6V0A1 transfected HeLa cells

1. Green - Atp6v0a1 antibody
2. Red- GM130 antibody (Golgi marker)
3. Merge with DAPI (blue)



For full product information, images and publications, please visit our [website](https://www.genetex.com).